

|               |                    |                            |                                |                            |                          |                                         |                        |                                  |                      |                                  |         |                                  |  |                                  |  |
|---------------|--------------------|----------------------------|--------------------------------|----------------------------|--------------------------|-----------------------------------------|------------------------|----------------------------------|----------------------|----------------------------------|---------|----------------------------------|--|----------------------------------|--|
| ALFA<br>INTES | Codice<br>FG519-01 |                            | Prodotto<br>Zerodue 10 ml      |                            | Formulazione<br>Collirio |                                         | Elaborato<br>Foglietto |                                  | Destinazione<br>Iraq |                                  | Braille |                                  |  |                                  |  |
|               | Carta              |                            | Dimensioni LxH<br>130 x 180 mm |                            | N. Col.<br>1             | Col. 1<br><div><div></div></div> P. 287 |                        | Col. 2<br><div><div></div></div> |                      | Col. 3<br><div><div></div></div> |         | Col. 4<br><div><div></div></div> |  | Col. 5<br><div><div></div></div> |  |
|               | Versione           | <div><div></div> 0.0</div> | <div><div></div> 0.1</div>     | <div><div></div> 0.2</div> |                          | <div><div></div> 0.3</div>              |                        | <div><div></div> 0.4</div>       |                      | <div><div></div> 0.5</div>       |         | <div><div></div> 0.6</div>       |  | <div><div></div> 0.7</div>       |  |
|               | Data               | 24.11.2025                 |                                |                            |                          |                                         |                        |                                  |                      |                                  |         |                                  |  |                                  |  |

## ENGLISH

# Zerodue

Ophthalmic Solution, 10 ml bottle

### Composition

Ophthalmic Solution, 10 ml bottle

### Composition

100 ml of solution contain:

Sodium hyaluronate 0.2 g, sodium chloride, magnesium chloride, potassium chloride, sodium phosphate monobasic, sodium phosphate dibasic, trisodium citrate, thimerosal 0.005 g, water for injection.

### Package

One box of Zerodue comprises 1 x 10 ml bottle.

### What Zerodue is

Zerodue is an ophthalmic solution.

### Why it should be used

Zerodue is an isotonic solution that, due to its sodium hyaluronate content (0.2%) and the presence of salts (magnesium, potassium, etc.) hydrates and lubricates the surface of the eye and reduce the risk of instability of the tear film. As it does not have a pharmacological effect, it protects the corneal epithelium from possible damage caused by specific conditions, such as:

- adverse weather and environmental conditions such as wind, cold, sun, salt water, air conditioning, smoke, smog and fine dust;
- specific ophthalmic stress conditions due to the use of contact lenses or prolonged work at video terminals or microscopes.

The lubricating action is also useful when, in specific physical conditions or after taking certain medications, tears are reduced, such as during the menopause and during treatment for diabetes.

### When it should not be used

Hypersensitivity to the constituents or other closely-related chemical substances.

### Precautions for use

Zerodue is for external use only. Zerodue must not

be used immediately before taking medications for treatment purposes or within 5 minutes of applying any pharmacological therapy via topical route.

Do not use the preparation whilst wearing contact lenses.

### How often and how long Zerodue should be used?

Zerodue should be dosed individually according to your condition and based on your eye doctor's recommendations. For general use, instil 1-2 drops of the solution into the conjunctival sac of each eye when needed, even several times a day.

### How to apply Zerodue

Instil one or two drops into each eye, allowing it to drop from above and preventing the bottle from touching your eyelid or eye surface.

### Side effects

In rare cases, a slight blurring of vision has been reported on instillation, which is due to the viscosity of the solution. This feeling, however, disappears very quickly. It is possible that, in rare cases, intolerance may occur (feeling of burning and/or irritation). Persistence of symptoms may be caused by the presence of the preservative. In this case, you should contact your eye doctor. Following the instructions contained in the information leaflet reduces the risk of side effects.

### Expiry and storage

Zerodue must be stored at a temperature of between 2°C and 25°C. The solution must be used within 4 weeks of opening the bottle. After this time, dispose of any remaining solution.

**Warning:** do not use the preparation after the expiry date stated on the package. Keep out of the reach of children.

CE 0373

STERILE A



### Manufacturer:

ALFA INTES  
Industria Terapeutica Splendore S.r.l.  
Via Fratelli Bandiera, 26  
80026 Casoria (Italy)

ALFA INTES