

Ocusilver

Preservative-free sterile gauzes

14 disposable sterile gauzes for topical periocular treatment in adults and children

INDICATIONS:

Ocusilver is indicated for eyelid and eyelash treatment both in adults and children, also in the presence of secretions resulting from inflammatory and/or infectious phenomena, such as blepharitis and blepharoconjunctivitis of varying nature (also in the presence of skin desquamation).

Its use is also particularly valid in the pre and/or post-operative phases of ophthalmic surgery or, in any case, in events resulting from trauma (wounds, abrasions, etc.). Ocusilver is also recommended for daily treatment of the periocular area to remove residues of cosmetic products from the eyelids and eyelashes in people who are predisposed to the onset of blepharitis and blepharoconjunctivitis of varying nature.

Owing to its particular composition, Ocusilver is also suitable for use by contact lens wearers.

PROPERTIES:

Sodium hyaluronate is a key component of the extracellular matrix of most tissues. Thanks to its ability to retain water and its film-forming power, it performs a moisturising and protective action.

All components of Ocusilver are useful for obtaining optimal protective and emollient action (chamomile and euphrasia), in full respect of the physiology of the particularly delicate and sensitive periocular area.

The product provides optimum chemical and physical compatibility with the tear film and it is able to remove secretions, any foreign bodies, dust, etc.

The distinctive composition of Ocusilver is indicated as an adjuvant in the treatment of inflammatory and infectious phenomena of varying nature in the periocular area.

COMPOSITION:

100% cotton gauze, soaked in 2.8 ml of the following solution: colloidal silver, sodium hyaluronate, chamomilla recutita, euphrasia officinalis, propylene glycol, phenethyl alcohol, PEG-40 castor oil, PEG-40 hydrogenated castor oil, sodium cocoyl wheat amino acids, ammonium glycyrrhizate, sodium bicarbonate, sodium chloride, aroma, polysorbate 20, ethylhexylglycerin, water.

INSTRUCTIONS FOR USE:

After washing your hands thoroughly, open the sachet and remove the gauze. Apply the gauze to the eyelid area, wiping it gently in order to adequately remove eye secretions and any other kinds of deposit (smog, fine dust and dust, etc.).

In case of particularly persistent secretions or occlusion of the lacrimal canaliculus, it is advisable to heat up the sealed packet for a few minutes in warm water (approximately 40°C) and then use the gauze as a warm compress, gently massaging the area.

In case of phlogosis, on the other hand, it is advisable to cool down the sealed packet in the refrigerator for a few hours and then apply the open gauze to the eye, keeping it in place for a few minutes as a cool compress.

WARNINGS:

Sterility of the product is only guaranteed if the packet is intact.

In case of infectious conditions, do not use the same gauze for both eyes. The gauze is disposable; reuse may be a vehicle for infection.

Store in a cool, dry place.

Keep out of the reach of children.

External use.

Discontinue treatment in case of evident irritative phenomena or hypersensitivity to the product.

Disposable.

Carefully read the instructions for use.



ALFA INTES

Industria Terapeutica Splendore S.r.l.

Via Fratelli Bandiera 26 - 80026 Casoria (NA)
Italy

CE 0425

STERILE R



ALFA INTES